

**İSTANBUL TEKNİK ÜNİVERSİTESİ**

**FEN-EDEBİYAT FAKÜLTESİ**

**GRADUATION PROJECT**



**MATLAB-based Nonlinear Optical Z-Scan Data Analysis Programming**

**Prepared By**

**MUHAMMED BEKTAŞOĞLU**

**090180132**

*Department : Physics Engineering Department*

**Advisor : Prof. Dr. Murat ERDEM**

**2023-2024 Spring Semester**

# CONTENTS

ABSTRACT	3
ACKNOWLEDGMENTS	4
FIGURES	5
ABBREVIATIONS AND SYBOMLS	6
1. INTRODUCTION	7
2. THEORETICAL PERSPECTIVE	8
2.1. Open Aperture Z-Scan	8
2.2. Closed Aperture Z-Scan	8
2.3. Two Photon Absorption	9
2.4. Three Photon Absorption	10
3. PROGRAM MANUAL	11
3.1. Installation	11
3.2. User Interface	15
3.3. Data Import	17
3.4. Data Analysis and Processing	20
3.5. Determination of Fitting Method and Parameters	22
3.6. Obtaining Graphical Output	24
3.7. Comprehensive Analysis and Export Features	26
4. ANALYSIS AND RESULTS	27
4.1. Application Scenario	27
5. CONCLUSION	30
REFERENCES	31
BIOGRAPHY	32

# ABSTRACT

This thesis presents a MATLAB-based application designed for analyzing experimental data from Z-scan measurements. Z-scan measurements are a critical technique for exploring nonlinear optical properties in materials, such as doped thin films and various types of glasses. The software's primary purpose is to aid in the importation and manipulation of raw data, allowing researchers to gain valuable insights into the optical characteristics of the materials under examination. The application features an intuitive user interface that simplifies the process of data manipulation, converting complex raw measurements into an easily analyzable format.

The software has advanced visualization tools that enable the graphical representation of scattering profiles. This allows for a direct visual assessment of sample behaviors under varying conditions. The application utilizes advanced fitting algorithms calibrated to align with theoretical models for phenomena such as open aperture Z-scan, closed aperture Z-scan, two-photon absorption and three-photon absorption processes. These algorithms accurately calculate critical material parameters, including the nonlinear absorption coefficient ( $\beta$ ) and the nonlinear refractive index ( $n_2$ ).

A MATLAB tool is available for this purpose, which integrates data importation, manipulation, and analysis in a user-friendly interface. This tool optimizes the research workflow. By streamlining these processes, the application not only enhances the efficiency and precision of nonlinear optical property studies but also contributes significantly to advancing the field of photonics research. The implementation of this application represents a significant step forward in the quest for a deeper understanding of nonlinear optical phenomena, offering researchers a powerful tool to expedite their discoveries in the realm of optical materials science.

# ACKNOWLEDGMENTS

I would like to express my deepest gratitude to Prof. Dr. Murat ERDEM, my thesis advisor, for his unwavering support and guidance throughout this research.

I cannot forget the incredible emotional support provided by my family, especially my grandparents, who have always believed in me.

MUHAMMED BEKTAŞOĞLU

July, 2024

# FIGURES

FIGURE 1: GENERAL INFORMATION .....	10
FIGURE 2: INSTALLATION DIRECTORY SELECTION .....	10
FIGURE 3: MATLAB RUNTIME INSTALLATION .....	11
FIGURE 4: USER AGREEMENT .....	11
FIGURE 5: INSTALLATION CONFIRMATION .....	12
FIGURE 6: INSTALLATION PROCESS .....	12
FIGURE 7: COMPLETION OF INSTALLATION .....	13
FIGURE 8: LAUNCHING THE APPLICATION .....	13
FIGURE 9: SAMPLE EXCEL SPREADSHEET FORMAT .....	14
FIGURE 10: THE EXCEL IMPORT BUTTON .....	15
FIGURE 11: FILE SELECTION DIALOG .....	15
FIGURE 12: IMPORTING DATA SET TO ORIGINAL TABLE .....	16
FIGURE 13: SCATTER PLOT OF IMPORTED DATA .....	16
FIGURE 14: DATA MANIPULATION CONTROLS – THE INTERFACE DISPLAYS THE MAIN CONTROLS FOR DATA MANIPULATION .....	17
FIGURE 15: EDITED DATA SET PREVIEW – THIS FIGURE SHOWS A SAMPLE OF THE EDITED DATA SET WITH TWO COLUMNS: Z (MM) AND INTENSITY .....	18
FIGURE 16: SELECTION MENU FOR FITTING METHODS .....	19
FIGURE 17: FITTING PARAMETERS INPUT – A GRAPHICAL USER INTERFACE PRESENTING A SET OF INPUT FIELDS FOR PARAMETERS USED IN FITTING PROCEDURES .....	20
FIGURE 18: SCATTERING AND FITTING SELECTION – AN INTERFACE CHECKBOX ELEMENT .....	21
FIGURE 19: EDITED DATA GRAPHS – THE LEFT GRAPH DISPLAYS AN EMPTY PLOT TITLED "EDITED DATA", PREPARED FOR DATA VISUALIZATION. THE RIGHT GRAPH SHOWS INTENSITY DATA AS A FUNCTION OF Z (MM), WITH A FITTED CURVE OVERLAID ON THE SCATTER PLOT, DEMONSTRATING THE RESULTS .....	21
FIGURE 20: COEFFICIENT OF DETERMINATION DISPLAY .....	22
FIGURE 21: CHART'S CONTEXT MENU .....	22
FIGURE 22: THE BUTTON IS USED TO EXPORT ALL PROCESSED DATA, INCLUDING THE ORIGINAL DATASET, THE MANIPULATED ENTRIES, AND THE RESULTS OF THE FITTING ANALYSIS, TO AN EXCEL FILE FOR FURTHER USE AND SHARING. ....	22
FIGURE 23: INTERFACE WITH "DATA EXPORT" AND "FIGURE EXPORT" BUTTONS, SELECTABLE DATA INPUTS, AND AN "OVERALL GRAPH" DISPLAYING TWO FIT LINES AND CORRESPONDING SCATTER DATA FOR COMPREHENSIVE DATA ANALYSIS AND VISUALIZATION. ....	23
FIGURE 24: EXCEL SPREADSHEET ILLUSTRATING THE EXPORT RESULTS FOR A DATA ANALYSIS PROCESS. ....	23
FIGURE 25: IMPORTING DATA FROM EXCEL AND ORIGINAL DATA SET .....	24
FIGURE 26: SCATTERING OF ORIGINAL DATA .....	24
FIGURE 27: DATA MANIPULATION FIELDS AND BUTTONS .....	24
FIGURE 28: CHOOSING THE METHOD BEST SUITED FOR THE DATA ANALYSIS NEEDS. ....	25
FIGURE 29: FIELDS TO ENTER VARIOUS PARAMETERS AND MENU ALLOWS USERS TO SELECT THE UNITS OF THE PARAMETERS .....	25
FIGURE 30: TABLE LISTS Z-VALUES AND THEIR CORRESPONDING INTENSITY VALUES, THE RESULTS FROM A FITTING ANALYSIS .....	25
FIGURE 31: THE GRAPH OF ANALYSIS, AND THE R SQUARE INDICATING THE FIT QUALITY OF THE MODEL TO THE DATA .....	25
FIGURE 32: REPRESENTING ANALYSES FROM MULTIPLE INPUT DATASETS SELECTED THROUGH CHECKBOXES UNDER 'INPUT TABS 1-5' ON THE LEFT SIDE OF THE INTERFACE. ....	26
FIGURE 33: PRESENTING PARAMETER VALUES OBTAINED FROM A DATA FITTING PROCESS. VARIOUS INPUT FIELDS AND RESULTS ARE CLEARLY DISPLAYED FOR USER ANALYSIS. ....	26

# ABBREVIATIONS AND SYMBOLS

$\beta$	:	Nonlinear absorption coefficient
$n_2$	:	Nonlinear refractive index
$Z$	:	Rayleigh Range
$L_{eff}$	:	Effective Length of the Sample
$\alpha$	:	Linear absorption coefficient
$L$	:	Thickness of the sample
$i$	:	Intensity of the incident laser beam
$\omega_0$	:	Beam waist
$\lambda$	:	Wavelength of laser

# 1. INTRODUCTION

Nonlinear optics is the study of the behavior of light in nonlinear media, where the material's response to light is not directly proportional to the light intensity. This field has revolutionized the way we manipulate light, leading to advancements in a wide range of applications from high-speed communication to medical imaging. Understanding the nonlinear optical properties of materials is crucial for the development of devices such as lasers, optical switches, and modulators.

## **Z-Scan Technique: A Tool for Nonlinear Optics Exploration**

The Z-scan technique is a simple yet powerful experimental method to measure the nonlinear optical properties of materials. It involves moving a sample through the focus of a laser beam and measuring the transmitted light intensity as a function of the sample position relative to the focal plane. This technique can distinguish between different types of nonlinearities, making it invaluable for characterizing materials for photonics applications.

The open aperture Z-scan method is used to measure nonlinear absorption, including phenomena like two-photon absorption and saturation absorption. In an open aperture setup, the entire transmitted beam is detected, allowing for the analysis of materials' absorptive nonlinearities without the influence of nonlinear refraction.

Contrastingly, the closed aperture Z-scan technique is utilized to assess nonlinear refraction. A small aperture is placed in the far field to partially block the transmitted beam. This setup enables the measurement of the phase distortion induced by the sample, providing insights into the refractive nonlinearities of the material.

Two-photon and three-photon absorptions are nonlinear optical processes where photons are absorbed simultaneously, transitioning an electron from a lower to a higher energy state. These phenomena are crucial for applications requiring controlled light-matter interactions, such as in optical limiting and photodynamic therapy.

This thesis introduces a MATLAB-based application designed to simplify the analysis of experimental data from Z-scan measurements, enhancing the study of nonlinear optical properties in materials. By offering an intuitive platform for data manipulation, and visualization, and fitting to theoretical models, this work contributes significantly to the field of photonics research. The application not only streamlines the experimental workflow but also provides a robust tool for the accurate characterization of nonlinear optical phenomena, thereby facilitating the advancement of optical materials science.

## 2. THEORETICAL PERSPECTIVE

### 2.1. Open Aperture Z-Scan

The Open-Aperture Z-Scan technique is used to measure the nonlinear absorption properties of materials by keeping the aperture fully open to capture all transmitted light. It is particularly effective at the focal point where light concentration is highest. Significant reductions in light transmission at this point indicate strong nonlinear absorption, such as two-photon absorption. This method is valuable for exploring the optical limiting capabilities of rare-earth ion-doped glasses and thin films, crucial for photonic and optoelectronic applications. Understanding these properties helps in designing optical limiters to protect sensitive equipment from intense laser pulses and in developing more efficient optical amplifiers. The open-aperture Z-Scan is noted for its sensitivity and versatility, making it a superior and essential tool in the study of optical nonlinearities compared to traditional methods. [1]

$$Y = 1 - 2 \left( \frac{\left(\frac{x}{z}\right)^2 + 3}{\left(\frac{x}{z}\right)^2 + 9} \right) \left( \left(\frac{x}{z}\right)^2 + 1 \right) \left( \frac{\beta I}{2} \right) L_{eff} \quad (2.1)$$

$$L_{eff} = (1 - e^{-\alpha L})/\alpha \quad (2.2)$$

$$z = \frac{\pi w_0^2}{\lambda} \quad (2.3)$$

### 2.2. Closed Aperture Z-Scan

The Closed-Aperture Z-Scan method is a precise experimental technique for determining the nonlinear refractive index ( $n_2$ ) of materials. In this method, a sample is positioned at the focal point of a laser beam that passes through an aperture. As the beam travels through the sample, it either expands or contracts due to changes in the refractive index, which are then detected at a second aperture set at twice the focal length from the lens. The technique utilizes changes in beam size to measure the nonlinear refractive index as the sample moves along the z-axis. Introduced by Sheik-Bahae et al. in 1990, it's known for its simplicity and effectiveness in quantifying nonlinear optical properties. The method operates on the principles of self-focusing and self-defocusing, which indicate whether a material has a positive or negative  $n_2$ . The detector measures changes in the transmitted beam's intensity to create a characteristic transmittance curve, which reflects the material's nonlinear properties. The curve's peaks and troughs help determine the sign and magnitude of  $n_2$ , providing valuable insights into the material's optical behavior. [2]

$$Y = 1 + \frac{4 \left( \frac{x}{z} \right) \left( \frac{2\pi}{\lambda} \right) n_2 L_{eff}}{\left( \left( \frac{x}{z} \right)^2 + 9 \right) \left( \left( \frac{x}{z} \right)^2 + 1 \right)} \quad (2.4)$$

$$L_{eff} = (1 - e^{-\alpha L})/\alpha \quad (2.5)$$

$$z = \frac{\pi w_0^2}{\lambda} \quad (2.6)$$



## 2.3. Two Photon Absorption

Two-photon absorption (TPA) is a distinct nonlinear optical phenomenon that involves a direct transition from the ground state to higher energy states facilitated by the simultaneous absorption of two photons. This process allows molecules to absorb energy from two photons concurrently, which is crucial for applications that require high precision and control over the absorption process.

In TPA, the typical energy states involved are quite different from those in single-photon absorption, making it a useful technique for studying materials with unique electronic properties. The challenge in distinguishing between two successive single-photon excitations and direct two-photon absorption lies in their similar end results; however, the mechanisms and implications of these processes can significantly influence the efficiency and resolution of optical systems. This makes understanding the nuances of TPA essential for advancing optical technologies. [3]

$$\frac{1}{(1+\beta L_{eff}(\frac{i}{(1+(\frac{x}{z})^2})))} \quad (2.7)$$

$$L_{eff} = (1 - e^{-\alpha L})/\alpha \quad (2.8)$$

$$Z = \frac{\pi w_0^2}{\lambda} \quad (2.9)$$

## 2.4. Three Photon Absorption

Three-photon absorption (3PA) is an advanced nonlinear optical phenomenon where three photons are simultaneously absorbed by a molecule, resulting in a transition from the ground state to an excited state. This process requires the simultaneous interaction of three photons whose combined energy is sufficient to reach an excited state without the need for intermediate state stabilization. 3PA is particularly significant in materials science and photonics, as it allows for deeper penetration and less scattering in biological tissues and other materials, making it ideal for deep-tissue imaging and photodynamic therapy.


$$\frac{1}{(1+2\beta L_{eff}(\frac{i}{(1+(\frac{x}{z})^2}))^2)} \quad (2.10)$$

$$L_{eff} = (1 - e^{-2\alpha L})/2\alpha \quad (2.11)$$

$$Z = \frac{\pi w_0^2}{\lambda} \quad (2.12)$$


## 3. PROGRAM MANUAL

### 3.1. Installation



*Figure 1: General Information*

1. Review the general information of the application and proceed by clicking the 'Next' button.



*Figure 2: Installation Directory Selection*

2. Select the directory where the installation will take place. It is recommended to leave it in the default **Program Files**.




Figure 3: MATLAB Runtime Installation

3. The MATLAB Runtime, necessary for running applications written in MATLAB, will be installed automatically. If MATLAB Runtime is already present on your system, this step will be skipped.






Figure 4: User Agreement

4. Accept the user agreement and click the 'Next' button.




*Figure 5: Installation Confirmation*

5. Confirm your selections and start the installation by clicking 'Begin Install'.




*Figure 6: Installation Process*

The installation will take a few minutes.



*Figure 7: Completion of Installation*

7. Once the installation is complete, exit by clicking the 'Close' button.



*Figure 8: Launching the Application*

8. To start the application, type "NonlinearOpticalZScanDataAnalysis" into the Windows search bar and click on it to launch.

## 3.2. Data Import

The *Data Import* functionality is a critical component of the program that allows for seamless integration of measurement data into the application. To ensure efficient data handling, the application supports importing data through the following steps:

### Multiple Independent Tabs:

The application is designed with five independent Input tabs, each capable of handling separate measurement data. This design enables concurrent processing of up to five distinct measurements within the same session, significantly enhancing productivity and multitasking capabilities.

### Exclusive Use of Excel:

Data importation relies exclusively on Microsoft Excel files. Users are required to prepare their measurement data in an Excel file format before initiating the import process.

### Data Format Specifications:


The data must be arranged in two columns, reflecting the dependency of Intensity on Z (mm) observed in Z-Scan configurations. It is crucial for the first column of data to reside in column A (representing Z(mm)), and the second column to be in column B (representing Intensity) within the Excel sheet.

### Headers Flexibility:

The presence or absence of axis headers in the Excel file is accommodated by the application. The import process is designed to recognize and correctly handle data with or without these headers, providing flexibility in data formatting.

### Single Observation Requirement:

Each Excel file should contain only the data for a single observation. This ensures clarity and prevents any potential confusion during the data analysis process. To maintain data integrity and simplicity, the Excel file should consist of a single worksheet. Data contained in secondary or multiple worksheets will not be recognized and hence should be avoided.



	A	B	C
1	x	y	
2	-75,2	96,958	
3	-74,575	98,099	
4	-73,95	97,719	
5	-73,325	98,099	
6	-72,7	98,479	
7	-72,075	97,719	
8	-71,45	97,719	

26	27	28	29
-60,2	96,578		
-59,575	97,719		
-58,95	98,099		
-58,325	98,859		


Figure 9: Sample Excel Spreadsheet Format

By adhering to the above guidelines, users can ensure a smooth and error-free data import experience into the application. To initiate the import, simply click on the green "Excel Import" button, conveniently located within each Input tab.



*Figure 10: The Excel Import Button*


Once the "Excel Import" button is activated, a file selector will appear, guiding users to select and import the desired Excel file from the computer.



*Figure 11: File Selection Dialog*

When the Excel file is selected, the data it contains is automatically transferred to the 'Original Data Set' section within the application. The import process seamlessly populates the table with data from the file and organizes them under predefined column headings such as 'Z (mm)' for distance measurement and 'Intensity' for observed values.

It is also instantly displayed as a scatter plot on the graph. This visual representation provides an immediate and intuitive view of the data points, allowing users to assess the distribution and key patterns in their dataset before proceeding with further analysis



16



### 3.3. Data Analysis and Processing

Following the successful import of data from Excel into the application, users can begin to manipulate and process the data for detailed analysis. This section of the manual guides you through the data manipulation features available after import.

#### Adjusting Data Alignment with Shift on x-direction:

The Z-Scan setup typically alters the distance starting from  $x = 0$  mm. If the central point of the graph is to be aligned at  $x = 0$  mm, users can achieve this by inputting a non-zero value in the "Shift on x-direction" field. This action will shift the  $Z(\text{mm})$  values in the first column of the table, effectively re-centering the dataset around the desired midpoint.

#### Data Truncation for Focused Analysis:

At times, observational data may exhibit undesired behavior beyond a certain range. To address this, users can specify their range of interest by replacing the default  $-\text{Inf}$  and  $\text{Inf}$  values in the "Data Between  $-\text{Inf} < z < \text{Inf}$ " section with the desired numerical limits. Only the data within the defined range will be retained for subsequent analysis, allowing users to exclude any portion of data that might distort the fitting process.

#### Normalization of Intensity Data:

By engaging the Normalized button, the Intensity data is normalized to scale between 0 and 1. This feature is particularly recommended if the data has not been previously normalized, ensuring uniformity and comparability within the dataset.




Figure 14: Data Manipulation Controls – The interface displays the main controls for data manipulation




### 3.4. Determination of Fitting Method and Parameters

Creating a precise fit for your data is crucial to accurately interpret your experimental results. In this section of the application, you will select a fitting method and define parameters that match your experimental setup.

#### Selecting the Fitting Method:

Firstly, choose the fitting method that aligns with the type of data you have. You'll see options like:

- Open Aperture Z-Scan
- Closed Aperture Z-Scan
- Two Photon Absorption
- Three Photon Absorption



*Figure 16: Selection menu for fitting methods*

After choosing the fitting method that best represents your experiment, the next critical step is to input the parameters that will be used to tailor the fit to your data.

#### Entering the Parameters:


Inputting parameters in the application is flexible, you do not need to fill in every field. The program is designed to work with either direct inputs from you or to calculate necessary values based on the information you provide. Here's how you decide which parameters to enter:

If you input the sample's thickness ( $L$ ) and the linear absorption coefficient ( $\alpha$ ), the application will automatically compute the effective length of the sample ( $L_{eff}$ ). In this case, you should not enter  $L_{eff}$  directly because the application calculates it for you.

Conversely, if you prefer to provide  $L_{eff}$  directly based on your calculations or knowledge, you should then leave the  $L$  and  $\alpha$  fields empty. The app is set up to prevent contradictory inputs.

For the wavelength ( $\lambda$ ) and beam waist ( $\omega_0$ ), inputting values will automatically restrict you from entering the Rayleigh range ( $Z$ ) as the program is designed to prevent the simultaneous input of these parameters. If you instead choose to enter a value for  $Z$  first, the application will then disable the entry fields for  $\omega_0$  to ensure the consistency of data required for the fitting process.

To perform a fit, the application must have an  $L_{eff}$  value. It can be either calculated by the app or input directly by you. Also, at least one of the following parameters must be provided for the fit to proceed:  $Z$  (Rayleigh Length),  $i$  (Intensity), or  $\beta$  (Nonlinear absorption coefficient). Without at least one of these, the fit cannot be performed.



The dialog box titled "Parameters for Fitting" contains two tabs: "Open Aperture Z-Scan" and "Closed Aperture Z-Scan". The "Open Aperture Z-Scan" tab is currently selected. It displays eight input fields arranged in two columns. Each field consists of a text label, a numeric input box, and a checkbox. The labels and their corresponding values are:  $L$  (0),  $\alpha$  (0),  $L_{eff}$  (0),  $i$  (0) in the left column; and  $\lambda$  (532),  $w_0$  (0),  $Z$  (0),  $\beta$  (0) in the right column. All checkboxes are currently unchecked.

Parameter	Value	Checkbox
$L$	0	<input type="checkbox"/>
$\alpha$	0	<input type="checkbox"/>
$L_{eff}$	0	<input type="checkbox"/>
$i$	0	<input type="checkbox"/>
$\lambda$	532	<input type="checkbox"/>
$w_0$	0	<input type="checkbox"/>
$Z$	0	<input type="checkbox"/>
$\beta$	0	<input type="checkbox"/>

Figure 17: Fitting Parameters Input – A graphical user interface presenting a set of input fields for parameters used in fitting procedures

### 3.5. Obtaining Graphical Output

After selecting the fitting method and entering the necessary parameters, it is crucial to ensure that data manipulation has been completed. This means that the "Edit Data Set" button must have been clicked to prepare the data in the "Edited Data Set" section.

If there are no issues with adhering to the rules for fit parameters, clicking on the **Get the scattering and fitting** checkbox initiates the application's calculation process.



Figure 18: Scattering and Fitting Selection – An interface checkbox element

The algorithm then proceeds to plot the optimal fit and scattering graph within the "Edited Data" graph and displays the results under the "Fit Data Set" section.

With these steps, the results are successfully obtained, showcasing the application's capability to seamlessly integrate data manipulation, parameter setting, and fitting processes to produce and display the fitting and scattering data comprehensively. This streamlined approach facilitates an intuitive and efficient analysis workflow, enabling users to derive meaningful insights from their data with ease.






Figure 19: Edited Data Graphs – The left graph displays an empty plot titled "Edited Data", prepared for data visualization. The right graph shows intensity data as a function of Z (mm), with a fitted curve overlaid on the scatter plot, demonstrating the results

Once the fitting process is complete, the coefficient of determination,  $R^2$ , which illustrates the compatibility between the scattering and fitting data, can be found in the lower right corner of the application. This  $R^2$  value is updated upon the completion of the fitting process, providing a quantitative measure of how well the fitted curve matches the original data points.




*Figure 20: Coefficient of Determination Display*

The context menu may be accessed by right-clicking on the chart. Within this menu, adjustments to the axis titles and chart title can be made. This context menu serves as a portal to various properties of the chart, typically initiated by a right-click action. It provides a convenient means to modify



*Figure 21: Chart's Context Menu*

After completing the data import and manipulation process, and once the fitting curve has been applied to your dataset, you will have a comprehensive set of data ready to be utilized for further analysis, reporting, or presentation. To facilitate this transition, the application provides an "Export" feature,



*Figure 22: The button is used to export all processed data, including the original dataset, the manipulated entries, and the results of the fitting analysis, to an Excel file for further use and sharing.*



## 4. ANALYSIS AND RESULTS

### 4.1. Application Scenario

In order for the analyses to be conducted, the Z(mm) and Intensity dataset obtained from the Z-Scan optical setup must be converted into a format that is acceptable by the application. To import Excel files in accordance with the structure that is accepted by the application, the **"Excel Import"** button is used. Refer to section [3.2 Data Import](#) to examine the suitability of the data structure.




Figure 25: Importing data from excel and original data set




Figure 26: Scattering of original data

Once the dataset has been incorporated into the application, data manipulation commences. Initially, the **'Shift on x-direction'** value is entered to adjust the positioning of the original data. To normalize the dataset and exclude certain data points, the endpoints of the data can be trimmed. To perform the editing process, the **'Editing'** button is finally pressed. For a comprehensive overview of data manipulation, please refer to section [3.3 Data Analysis and Processing](#).




Figure 27: Data Manipulation fields and buttons

Prior to commencing the analytical process, it is necessary to select the most appropriate fitting method. Once this has been determined, the parameters of the system must be entered. In order for the analysis to be performed, it is essential to input at least the  $L_{eff}$  value (or L and alpha,



if calculable). The program has been designed to calculate  $Z$  (Rayleigh Length),  $i$  (Intensity), or  $\beta$  (Nonlinear absorption coefficient) values; however, at least one of these parameters must be entered. If the intensity or Rayleigh length is known, entering the  $i$  (Intensity) value, or  $Z$  value, will facilitate the calculation of the remaining parameter.

However, if the application attempts to calculate two unknown values, the processing speed will slow down due to the principle of nested loops. Calculating a single unknown is much simpler. The calculation of parameters is mandatory, and there may be different calculation scenarios. For detailed information about the fitting method and parameters, it is recommended that the reader review section [3.4. Determination of Fitting Method and Parameters](#).




Figure 28: Choosing the method best suited for the data analysis needs.




Figure 29: Fields to enter various parameters and menu allows users to select the units of the parameters

Once the requisite fit parameters and units have been entered and the **"Fitting"** button has been pressed, the application initiates a series of processes, which run in the background depending on the chosen method and parameters. Once the application completes the analysis, a fit graph is generated, and the data from the line on the graph is presented to the user as the **"Fit Data Set"**

To view the entire data set together, the **"Export"** button can be used to download the data.




Figure 31: The graph of analysis, and the R square indicating the fit quality of the model to the data.

Z (mm)	Intensity	$\pm$
-59.5750	0.9727	
-58.9500	0.9722	
-58.3250	0.9717	
-57.7000	0.9711	
-57.0750	0.9705	
-56.4500	0.9699	
-55.8250	0.9693	
-55.2000	0.9687	
-54.5750	0.9681	
-53.9500	0.9674	
-53.3250	0.9667	
-52.7000	0.9660	

Figure 30: Table lists Z-values and their corresponding intensity values, the results from a fitting analysis

After working on different datasets, switch to the **Outputs** tab to compare all results from Input Tabs 1-5. The desired datasets are selected from the **Input Tabs** and displayed overlaid on the **Overall Graph**.




Figure 32: Representing analyses from multiple input datasets selected through checkboxes under 'Input Tabs 1-5' on the left side of the interface.

Additionally, the results of the obtained parameters can be accessed from the area located in the bottom right of the Outputs tab.

Input 1	Input 2	Input 3	Input 4	Input 5
Thickness of the material (L)	0	Wavelength of laser ( $\lambda$ )	0	
Linear absorption coefficient ( $\alpha$ )	0	Beam Waist ( $w_0$ )	0	
Effective Length of the Sample ( $L_{eff}$ )	0.0039	Rayleigh length (Z)	7.6	
Intensity of the incident laser beam (i)	1e+07	Nonlinear absorption coefficient ( $\beta$ )	4.708e-05	
D	0	Reduced Chi-Square	0.002704	
Units of the parameters :	millimeter    Watt centimeter    Watt meter    Watt	R-Square (COD)	0.9128	
		Adjusted R-Square	0.9122	

Figure 33: Presenting parameter values obtained from a data fitting process. Various input fields and results are clearly displayed for user analysis.

## 5. CONCLUSION

The development and implementation of the MATLAB-based application described in this thesis significantly enhances the analysis of Z-scan measurement data, a cornerstone technique in the study of nonlinear optical properties. The comprehensive suite of features, including data importation, advanced fitting algorithms, and intuitive visualization tools, enables the application to successfully address the complexities involved in decoding the nonlinear optical behavior of various materials.

The software offers a streamlined, user-friendly platform that not only simplifies the research process but also enables more accurate and efficient exploration of critical parameters such as the nonlinear absorption coefficient and the nonlinear refractive index. The inclusion of both open and closed aperture Z-scan techniques, in addition to the capabilities for analyzing multi-photon absorption phenomena, ensures that researchers can conduct a thorough investigation of material properties under diverse experimental conditions.

This thesis has demonstrated that the application considerably reduces the time and effort required to process and interpret experimental data, allowing researchers to focus more on the implications of their findings rather than the intricacies of data handling. Moreover, the graphical output generated by the software facilitates a better understanding and communication of results.

# REFERENCES

- [1] Theoretical-experimental method of nonlinear optics: the Z-scan technique South Florida Journal of Development, Miami, v.4, n.4. p. 1807-1827, 2023. ISSN 2675-5459
- [2] Sheik-Bahae, M., Said, A. A., & van Stryland, E. W. (1990). Sensitive measurement of optical nonlinearities with a single beam Z-scan. Optics Letters, 14(11), 955-957.
- [3] Atomically Precise Metal Nanoclusters. (2022). In section 4.6.1 Two-photon absorption and emission.
- [4] Femtosecond and nanosecond nonlinear optical properties of alkyl phthalocyanines studied using Z-scan technique



## BIOGRAPHY

Muhammed was born in Istanbul in the year 2000. He spent his childhood and teenage years and completed his primary, secondary, and high school education in Sarıyer, İstanbul. In 2018, he commenced his studies in Physics Engineering at Istanbul Technical University, where he had the opportunity to merge theoretical knowledge with practical application.

Additionally, he participated in the Erasmus exchange program at Umeå University in Sweden, in the third year of the study, further broadening his academic and cultural perspectives.

During his university studies, Muhammed gained professional experience by working as a full-stack web developer and system analyst at Turkish Airlines and Amadeus IT Service. These experiences helped enhance his passion and skills in technology and software.

Currently, in the final year of his undergraduate degree, he is engaged in significant research in the field of optics, collaborating with his advisor, Prof. Dr. Murat Erdem.